



**Price List and Order Form**

Effective April 1, 2007 Revised March 17, 2009

For ease of ordering, contact customer Service Toll Free at:

**\*\*Recent updates are in RED**

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

**Yes, you can.**

Company Name \_\_\_\_\_ Account # \_\_\_\_\_

Phone \_\_\_\_\_ P.O. # \_\_\_\_\_ Date \_\_\_\_\_

Ship to Name & Address \_\_\_\_\_

**Special Note**

**HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.**

**FRAME SELECTION**

**POWER CHAIR TYPE**

- M91R M91 w/SureStep w/MK6i™ Elect.(K0848) .....5,600.00  
(Push button style seat positioning strap included)
- M91-C M91 Complete w/ SureStep Base Only, with  
MK6i electronics .....5,255.00
- M91-M M91 Motion Concepts w/SureStep Base Only, with  
MK6i electronics .....5,255.00
- M91SEAT Pronto M91 Seat Only .....740.00

**TRANSPORT TIE DOWN**

- TRRO Transport Ready Option .....300.00
- TRBKTS Wheelchair Transport Brackets .....250.00

**\* Transport Tie Down**

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNAWC Vol 1. Section 19.5.3 Frontal Impact Test requirements for wheelchairs with a 168 lb. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

\*Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

**TRRO FEATURES FOR M91 WITH SURE STEP**

FEATURES	
<ul style="list-style-type: none"> <li>• Pronto M91R w/MK6i</li> <li>• 300 lb or 400 lb User Weight Capacity</li> <li>• Adjustable Rehab Seat with Adj Seat and Tilt Back Angle</li> <li>• Any Seat Width</li> <li>• Any Seat Depth</li> <li>• Sectional Rear Uprights</li> <li>• 18" – 20" Back Height</li> <li>• Any Arm Type</li> <li>• Black Nylon Upholstery</li> <li>• 14" x 3" Flat Free Tires</li> <li>• Any Shroud Color</li> <li>• Footboard</li> </ul>	<ul style="list-style-type: none"> <li>• MPJ+™</li> <li>• All PSR/PSF Joysticks</li> <li>• Left Handed Mounted Joystick Option</li> <li>• "T" Handle Flexible Joystick Ext.</li> <li>• Straight Handle Flexible Joystick Extension</li> <li>• MK6i Remote Programmer</li> <li>• Reinforced Back Upholstery</li> <li>• 22NF Batteries</li> <li>• Tire Pump</li> <li>• Lester 8 amp Battery Charger</li> </ul>

**SEAT SELECTION**

- ADJASBA Seat w/ Adj. Seat Tilt & Back Angle ..... **Standard**
- CTS Contoura™ Seating ..... 395.00 <sup>1</sup>
- OMCU Omit Contoura Seat Cushion ..... Less 30.00
- SEMIRB Van Seat - Semi Recline with Solid Base ..... 70.00

**MOTOR PACKAGES (USER WEIGHT LIMIT)**

- U300 User Weight Up to 300lbs (STDPKG) (K0848) **Standard** <sup>2</sup>
- U400 User Weight 301 – 400 lbs (HDPKG) (K0850) ... 400.00

**ADJUSTABLE ASBA SEAT**

**ADJUSTABLE SEAT WIDTH RANGE**

- 1620W 16" to 20" Seat Width Range ..... No Charge
- 2024W 20" to 24" Seat Width Range ..... No Charge

**SEAT WIDTH**

- W16 16" Seat Width ..... No Charge
- W17 17" Seat Width ..... No Charge
- W18 18" Seat Width ..... **Standard**
- W19 19" Seat Width ..... No Charge
- W20 20" Seat Width ..... No Charge
- W21 21" Seat Width ..... 329.00
- W22 22" Seat Width ..... 329.00
- W23 23" Seat Width ..... 329.00 <sup>3</sup>
- W24 24" Seat Width ..... 469.00 <sup>3</sup>

**ADJUSTABLE SEAT DEPTH RANGE**

- 1619D 16" to 19" Seat Depth Range ..... No Charge
- 1922D 19" to 22" Seat Depth Range ..... No Charge

**SEAT DEPTH**

- D16 Seat Depth 16" ..... **Standard**
- D17 Seat Depth 17" ..... No Charge
- D18 Seat Depth 18" ..... No Charge
- D19 Seat Depth 19" ..... No Charge
- D20 Seat Depth 20" ..... 529.00
- D21 Seat Depth 21" ..... 529.00
- D22 Seat Depth 22" ..... 529.00

**UPHOLSTERY COLORS**

- U240 Black Nylon Upholstery ..... No Charge

**FOOTNOTES**

<sup>1</sup> Selected Seat Depth on Contoura assumes 1" of back cushion compression. For Example: When choosing 18" depth the measurement from the front of the back cushion to the front of the seat pan equals 17". The contoura back cushion is 2" thick.

<sup>2</sup> Standard motor package includes 4-pole motors with a weight limit of 300lbs.

<sup>3</sup> Only available with U400 motor package

# Invacare Pronto M91 with SureStep (Rehab Version)

## Accessories and Options

### BACK SELECTIONS

#### SEAT BACK HEIGHT

- BH16 Back Height 16" ..... No Charge
- BH17 Back Height 17" ..... No Charge
- BH18 Back Height 18" ..... No Charge
- BH19 Back Height 19" ..... No Charge
- BH20 Back Height 20" ..... No Charge

#### CONTOURA

- TSBH25 Back Height 25" ..... No Charge
- TSBH295 Back Height 29.5" ..... No Charge

#### SEAT BACK TYPE

- SECT Sectional Rear Uprights ..... No Charge

### ARM SELECTIONS

#### ARM TYPE

- 28F Desk Length Fixed Height Conv. .... **Standard**
- 29F Desk Length Adj Height Conv (E0973) ..... 230.00
- 18F Full Length Fixed Height Conv. .... No Charge
- 19F Full Length Adj Height Conv (E0973) ..... 230.00

#### MIXED HEIGHT ADJUSTABLE ARM TYPE (EACH)

- 19FR Full Length Adj Hgt Right (E0973) ..... 115.00
- 29FR Desk Length Adj Hgt Right (E0973) ..... 115.00
- 19FL Full Length Adj Hgt Left (E0973) ..... 115.00
- 29FL Desk Length Adj Hgt Left (E0973) ..... 115.00

### VAN SEAT – WITH SOLID BASE

#### VAN SEAT

- VW16 16 Wide Van Seat ..... No Charge
  - VD16 16 Depth ..... No Charge
  - VD17 17 Depth ..... No Charge
  - VD18 18 Depth ..... No Charge
- VW18 18 Wide Van Seat ..... Standard
  - VD16 16 Depth ..... Standard
  - VD17 17 Depth ..... No Charge
  - VD18 18 Depth ..... No Charge
- VW20 20 Wide Van Seat ..... 329.00
  - VD18 18 Depth ..... No Charge
  - VD19 19 Depth ..... No Charge
  - VD20 20 Depth ..... 529.00
- VW22 22 Wide Van Seat ..... 529.00
  - VD18 18 Depth ..... No Charge
  - VD19 19 Depth ..... No Charge
  - VD20 20 Depth ..... 529.00

#### VAN UPHOLSTERY

- GV Grey Vinyl ..... No Charge

#### ARM TYPE

- VAH Full Length Adjustable Height (E0973) ..... 230.00
- VFH Full Length Fixed Height ..... Standard
- VAHD Desk Length Adjustable Height (E0973) ..... 230.00
- VFHD Desk Length Fixed Height ..... No Charge

### POSITIONING STRAPS & WHEEL OPTIONS

#### SEAT POSITIONING STRAPS

- 1515 Push Button Style Seat Positioning Strap ..... **Standard**
- 1311BK Airline Buckle Seat Positioning Strap ..... No Charge

#### CHEST POSITIONING STRAPS

- 7311BK Airline Buckle Chest Pos. Strap (E0960) ..... 100.00
- 9515 Push Button Style Chest Pos. Strap (E0960) ..... 100.00
- 7321BK Hook & Loop Chest Pos. Strap (E0960) ..... 100.00

### SHROUD COLORS

- BL Deep Blue Pearl ..... No Charge
- RD Candy Red Pearl ..... No Charge
- MB Magic Black ..... No Charge

### DRIVE WHEEL TYPE

- DW143 14" X 3" Pneumatic Tires ..... **Standard**
- DW143FF 14" X 3" Flat Free Tires ..... 149.00

### ACCESSORIES

#### BACK ACCESSORIES

- L8601 Omit Back Upholstery ..... No Charge
- RBP Reinforced Back Upholstery ..... No Charge

#### ARM ACCESSORIES

- L8401 Omit Arms ..... Less 20.00
- L8405 Omit Arm Pads ..... Less 10.00

#### MISCELLANEOUS ACCESSORIES

- SEAT02 Oxygen Holder (ASBA Seat Only) (E2208) ..... 150.00
- HFO2-ML6 HomeFill™ O2 Holder-ML6 Only (E2208) ..... 150.00
- 1209 Tire Pump ..... 35.00
- LESTR Lester 8 Amp Battery Charger (E2367) ..... 320.00
- VANO2 Oxygen Holder (Van Seat Only) (E2208) ..... 150.00
- CTE Single Mode Charger (E2366) ..... 205.00

### BATTERY SELECTION

#### BATTERIES

- 22NFINSTAL 22NF Battery (On-Chair)(1 pair)(E2361) ..... 300.00

### MK6I ELECTRONICS

#### ELECTRONIC CONTROLLERS

- MK690 MK6i 90 Amp Controller ..... **Standard**

#### STANDARD DRIVER CONTROL

- SPJM6 SPJ™+ Single Drive Joystick ..... **Standard**

#### BASE ONLY SPJ+ OPTIONS <sup>4</sup>

- SPJAP SPJ+ with Powered Seating Switch ..... **Standard**
- SPJAPM Single Actuator Operation thru SPJ+ (E2310).. 1,300.00 <sup>5</sup>

### FOOTNOTES

<sup>4</sup> Base Only Systems for Multiple Actuator Powered Seating Systems will require MPJ+, PSRM6, PSFM6, or MK6i Display and Alternative Driver Control.

<sup>5</sup> Allows Actuator Operation through Driver Control. For Single Actuator Powered Seating Systems.

# Invacare Pronto M91 with SureStep (Rehab Version)

## Accessories and Options

### EXPANDABLE CONTROLLER SYSTEMS

- PWH Harness req. for Expandable Electronics (E2313) 350.00<sup>6</sup>
- MPJM6U MPJ™+ Multiple Drive Joystick (E2377) ..... 750.00
- PSFM6 Personal Joystick, Inductive on Front ..... 1,130.00<sup>7</sup>
- PSRM6 Personal Joystick, Inductive on Rear ..... 1,130.00<sup>7</sup>
- PBOD Push Button On/Off w/PB Drive Select .... No Charge<sup>8</sup>
- PBSS Push Button On/Off w/Speed Select ..... No Charge<sup>8,9</sup>
- OR On/Off Mounted Right ..... No Charge<sup>8</sup>
- OL On/Off Mounted Left ..... No Charge<sup>8</sup>
- ASLRDYM6 MK6 Display Only (No Driver Control) (E2377) 750.00<sup>10</sup>

**FOR YOUR ASL DRIVE SYSTEM NEEDS  
PLEASE REFER TO THE FOLLOWING:  
ASL SINGLE INVOICE PRICE LIST FORM # 02-097**

- 1500M6 Proportional RIM Head Control (E2327) ..... 2,875.00
- SMHD Small Headrest ..... No Charge
- LGHD Large Headrest ..... No Charge
- 1558MM6 Compact Joystick w/Display (E2373) ..... 1,130.00

### ALTERNATE DIGITAL CONTROL SYSTEMS

- SNPM6 Sip-N-Puff Interface Box (E2325) ..... 1,500.00<sup>11</sup>
- 5018M6 Wafer Board & D/I Box (E2322) ..... 1,550.00<sup>11</sup>
- 5020M6 Mini Tash Joystick & D/I Box (E2321) ..... 1,750.00<sup>11</sup>

### SIP N PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit (E2326) .. 380.00

### ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount ..... 770.00<sup>12</sup>
- 5018 Wafer Board ..... 310.00<sup>12</sup>
- 5020 Mini Tash Joystick ..... 405.00<sup>12</sup>
- PACM6 Proportional Attendant Control ..... 850.00<sup>12</sup>
- 1552M Digital Attendant Control ..... 850.00

### JOYSTICK MOUNTING HARDWARE

- QL Quad Link Retract Joystick Mnt SPJ+ (E1028) .... 230.00<sup>13</sup>
- QLM6 Quad Link Retract Joystick Mnt MPJ+ (E1028) .. 279.00<sup>13</sup>
- ARM250 Stealth Height Adjustable Joystick Mount ..... 210.00
- RIMHW Cane Mounted Head Control Hardware (E1028) . 219.00
- GATMPJ6 Gatlin Midline for MPJ (E1028) ..... 600.00
- GAT1812 Gatlin Midline for Compact Joystick (E1028) .... 600.00
- PKG32669 Complete Bib Assembly ..... 265.00
- RCM Rear Cane Mount for Compact Joystick ..... 395.00
- LEFT Left Handed Mounted Joystick ..... No Charge

### FOOTNOTES

- <sup>6</sup> Harness is required to interface joystick or MK6i display with expandable controllers. Option is located under Electronic Accessories in configurator
- <sup>7</sup> TOD-SS: Toggle On/Off- Drive Select, Speed Potentiometer Standard.
- <sup>8</sup> PSRM6 Only.
- <sup>9</sup> Drive Select through mode switch, or mode switch and Joystick Command.
- <sup>10</sup> Display Only, No Driver Control. Allows choice from order form.
- <sup>11</sup> Must also use either a Multiple Drive Joystick or MK6i Display.
- <sup>12</sup> Requires SNPM6 Sip-N-Puff Interface Box.
- <sup>13</sup> QL and QLM6. ARM250 does not come standard with Quad Links on Pronto Series.

### ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer w/Pro Memory Card ..... 475.00
- 1813SD Pro Memory Card (USB Ready) ..... 90.00
- 1560 "T" Handle Flexible Joystick Ext. (E2323) ..... 77.00
- 1561 Straight Handle Flexible Joystick Ext. (E2323) ..... 77.00
- 1826 Chin Cup (E2324) ..... 48.00
- A24VPS 24V MK6 Auxiliary Power Source ..... 200.00
- AUX12M6 Auxiliary Module For 1 & 2 ..... 795.00
- AUX34M6 Auxiliary Module For 3 & 4 ..... 795.00
- MPY Monoport "Y" Cable ..... 15.00

### ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch – Black ..... 65.00
- ASL304A Light Wobble Switch w/611 Mount ..... 320.00
- ASL202J Fiber Optic Switch ..... 610.00
- ASL208J Adjustable Proximity Switch ..... 550.00
- ASL504M6 Remote Emergency Stop Switch ..... 1,185.00

### FRONT RIGGING OPTIONS

#### FRONT RIGGING

##### FOOTBOARD

- 90FB Footboard ..... **Standard**
- L90FB Omit Footboard ..... **Less 40.00**

##### FOOTREST/LEGREST (HEMI)

- PHW93 Swingaway Footrests 13½" – 19" ..... 165.00
- PH904A Swingaway ELR 16" – 19¾" (E0990) ..... 165.00
- PHAL4A Smartleg® Art. Legrests 17¼" – 21½" (K0053) .. 225.00

##### FOOTPLATE OPTIONS

- 1350 Extra Large Aluminum Footplates ..... 115.00
- 1450 Aluminum Footplates Deluxe ..... **No Charge**
- 1650 Aluminum Footplates ..... **No Charge**
- AT5543 Adjustable Angle Flip Up Footplates ..... 165.00

##### FOOTREST/LEGREST ACCESSORIES

- LPT2 Longer Pivot Slide Tube For 1350 ..... 83.00
- A1337 Calf Strap (K0038) ..... 27.00
- IMPACTG Impact Guards ..... 51.00
- ALPT Longer Pivot And Slide Tube ..... 84.00
- ALPT4 Longer Pivot And Slide Tube 4" ..... 84.00

### 5 STAR SERVICE PLAN

#### 5 STAR SERVICE PLAN

- 5STAR1 Invacare 5 Star Service Plan - Guarantees predictable service costs for one year on Invacare Consumer Power and Scooters ..... 150.00<sup>14</sup>

### SEATING & POSITIONING

**FOR YOUR SEATING NEEDS  
PLEASE REFER TO THE FOLLOWING**

**INFINITY® SEAT AND BACK CUSHIONS  
See Price List Form # 01-351**

**CONVENTIONAL PRODUCTS  
See Price List Form # 95-82**

**ABSOLUTE™, LAP TRAYS, AND PERSONAL BACK PRODCUTS  
See Price List Form # 01-207**

### FOOTNOTES

- <sup>14</sup> Visit [www.invacare.com/ipartner](http://www.invacare.com/ipartner) for more information.

**NOTE: All specifications and dimensions are approximate.  
For the most current pricing information go to [www.invacare.com](http://www.invacare.com).**

# Invacare Pronto M91 With SureStep (Rehab Version)

## Accessories and Options

### Memory Card – BASIC

- Standard with ALL Expandable Control Systems (MPJ+, PSR+, PSF+, MK6i Display).
- Used to Back Up/Restore Programmed Settings for Only One Chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not contain Advanced Diagnostics/Help Screens/File Structure.
- Not compatible with SPJ+ Joysticks.

### Memory Card – PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, Help Screens, File Storage/Retrieval.
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

## SPECIFICATIONS

<b>Dimensions</b>	
Lengths (with footboard folded)	39"
Lengths (with footboard down)	40" – 42"
Widths (excluding Joystick)	25.9"
Heights	37" – 40"
<b>Seat Dimensions</b>	
ADJASBA Seat	
Widths	16" – 24"
Depths	16" – 22"
Back Height	16" – 20"
<b>Seat-to-Floor Height</b>	
ADJASBA	18.5" – 21"
<b>Weight</b>	
With batteries	200 lbs.
Without batteries	150 lbs.
<b>Weight of Components Disassembled</b>	
Battery and Cover Assembly (2 per unit)	25 lbs.
Battery tray, charger, controller shrouds	21 lbs.
Side frame, motor gearbox, side shrouds (2 per unit)	39 lbs. (ea.)
Footboard and mounting tube	5 lbs.
<b>Wheel and Caster sizes</b>	
Front Casters	6" x 2"
Drive Wheels	14" x 3"
Rear Casters	6" x 2"
Maximum Speed	6 mph (Standard) 5 mph (Heavy duty)
Incline Capability	9°
Weight Capacity	300 lbs. (400 lb. option)
<b>Turning Radius</b>	
Front (with footboard down)	19.5"
Rear	21.5"
Ground Clearance	3.25"
<b>Warranties</b>	
Base frame	60 months
Seat Frame	36 months
Electronics	12 months
Motor/gearbox	12 months



**Yes, you can.**

### Invacare Corporation

#### USA

One Invacare Way  
Elyria, OH USA  
44035-4190  
800-333-6900

[www.invacare.com](http://www.invacare.com)

All rights reserved. Trademarks are identified by the symbols <sup>TM</sup> and <sup>®</sup>. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Products designated by ASL are Manufactured by Adaptive Switch Labs, Inc. Therafin is a registered trademark of W/R Medical Electronics Co. Specifications and prices are subject to change without notification.

©2007 Invacare Corporation

Form No. 04-211 Rev 0309